




Rayon Tipped Applicator

Polystyrene handle



Description

This specialty swab is made with USP grade rayon fiber that is carded on-site to maintain consistency. This swab is ideal for diagnostic specimen collection, wound treatment and application of medicines. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive.

Dimensions

3" Overall Length (76.2 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 2.937" Length (74.6 mm)

Qualities

- Medical grade quality
- Absorbent tip
- Soft smooth tip
- Semi-flexible polystyrene handle
- Made in USA

Applications

- All purpose general swab
- Application of topical dressings
- First aid use
- Applying and removing waxes in health/beauty applications
- Cleaning sensors and components

Ref #	Packaging
803-PR Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)



Puritan®