



# HydraFlock® Flocked Swab - Dry Transport System

Elongated tip, polystyrene handle with breakpoint at 100 mm



## Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kit. The patented HydraFlock® flocked swab provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

## Dimensions

### Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.00" (152.4 mm)

### Tip

- 0.219" Dia. (5.563 mm)
- 0.685" Length (17.399 mm)

### Handle

- 0.100" Dia. (2.54 mm) tapered to 0.067" (1.702 mm)
- 5.96" Length (151.384 mm)

## Qualities

- Medical grade quality
- Highly absorbent capacity
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

## Applications

- Rapid diagnostic testing
- Forensic evidence collection
- Specimen collection
- For buccal cell collection
- DNA testing

Ref #	Packaging
25-3706-H BT Sterile	5/100/1 (5 boxes, 100 tubes per box, 1 swab per tube)



**Puritan®**