



# Popule® Self-Saturating Foam Swab

Handle filled with 91% and 9% DI water, foam tip



-Self-saturating foam tip

## Description

This pre-filled popule is filled with 91% IPA and 9% distilled deionized water. It has a reticulated 100PPI foam tip which is bonded to a clear polypropylene handle with a hot melt adhesive. Snap the handle to release the liquid into tip - neat and clean for your convenience. Ready for use in the field for the collection of Trace DNA

## Dimensions

3.312" Overall Length (76.2 mm)

### Tip

- 0.438" Dia. (11.125 mm)
- 1.00" Length (25.4 mm)

### Handle

- 0.270" Dia. (6.858 mm)
- 3" Length (76.2 mm)

## Qualities

- Self-saturating foam tip
- Non-glass handle
- Handle remains dry and prepping solution stays off fingers
- No cross contamination
- Instructions for use included in each box
- Made in USA

## Applications

- Ideal for forensics/DNA
- Ready for use in the field for the collection of Trace DNA
- General removal of residues
- Critical environment applications

Ref #	Packaging
25-4660 Sterile	10/40 (10 boxes per case, 40 per box)



**Puritan®**