



Infant Tongue Depressor

Wood



- Biodegradable wood handle

Description

This infant size blade is made of Northern Maine white birch which is naturally biodegradable. It is available non-sterile or sterile depending on the application. Primarily used in oral patient care on infants and small children.

Dimensions

4.5" Overall Length (114.3 mm)
0.375" Width (9.525 mm)

Qualities

- Medical grade quality
- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges
- Made in USA

Applications

- Oral patient care and examinations
- Veterinary examinations
- Applying or mixing medicines

Ref #	Packaging
700 Non-sterile	10/1000 (10 boxes per case, 1,000 per box)
25-700 Sterile	10/100/1 (10 boxes of 100 individually wrapped packages, 1 p/ pkg)



Puritan[®]