

704 / 704 BULK

Standard Tongue Depressor

Non-sterile, standard size wood blade, available in bulk

DESCRIPTION

This standard size blade is made of Northern Maine white birch which is naturally biodegradable. Used in oral patient care and examinations, health and beauty applications, and a variety of industrial uses.

FEATURES & BENEFITS

- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges
- #704 – Puritan brand
- Made in USA

APPLICATIONS

- Oral patient care and examinations
- Veterinary examinations
- Health and beauty; waxing
- Stirring stick
- Applying epoxy or resins
- Industrial use

Item No.	Dimensions length x diameter	Packaging pack, weight
non-sterile sterile		
704	152.4 x 17.475 mm 6 x 0.688"	10/500,* 35 lbs
704 BULK		3,700 38 lbs (bulk)

*10 boxes per case, 1,000 per box

Products may not appear at actual size. See chart or our website for accurate measurements.

