




# Adult Tongue Depressor

Wood



- Biodegradable wood handle

## Description

This adult size tongue blade is made of Northern Maine white birch which is naturally biodegradable. Available sterile or non-sterile depending on the application. Primarily used in oral patient care and examinations.

## Dimensions

6" Overall Length (152.4 mm)  
0.75" Width (19.05 mm)

## Qualities

- Medical grade quality
- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges
- Made in USA

## Applications

- Oral patient care and examinations
- Veterinary examinations
- Applying or mixing medicines
- Health and beauty applications - waxing
- Stirring stick
- Applying epoxy or resins

Ref #	Packaging
705 Non-sterile	10/500 (10 boxes per case, 500 per box)
25-705 Sterile	10/100/1 (10 boxes of 100 individually wrapped-packages, 1 p/pkg)



**Puritan®**