

706

Standard Tongue Depressor

Non-sterile, standard size wood blade

DESCRIPTION

This standard size blade is made of Northern Maine white birch which is naturally biodegradable. Primarily used in oral patient care and examinations.

FEATURES & BENEFITS

- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges
- Made in USA

APPLICATIONS

- Oral patient care and examinations
- Veterinary examinations
- Health and beauty; waxing
- Stirring stick
- Applying epoxy or resins

Item No.	Dimensions length x diameter	Packaging pack, weight
non-sterile sterile		
706	152.4 x 17.475 mm 6 x 0.688"	12/100,* 9 lbs

*12 boxes per case, 100 per box

Products may not appear at actual size. See chart or our website for accurate measurements.

