

709

Standard Tongue Depressor

Non-sterile, standard size wood blade, individually wrapped

DESCRIPTION

This standard size tongue blade is made of Northern Maine white birch which is naturally biodegradable. Primarily used in oral patient care and examinations.

FEATURES & BENEFITS

- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges
- Individually wrapped
- Made in USA

APPLICATIONS

- Oral patient care and examinations
- Veterinary examinations
- Health and beauty; waxing
- Stirring stick
- Applying epoxy or resins

Item No.	Dimensions length x diameter	Packaging pack, weight
non-sterile sterile		
709	152.4 x 17.475 mm 6 x 0.688"	10/250,* 21 lbs

*10 boxes per case, 250 per box, individually wrapped

Products may not appear at actual size. See chart or our website for accurate measurements.

