



Standard Tongue Depressor

Wood



- Biodegradable wood handle

Description

This standard size tongue blade is made of Northern Maine white birch which is naturally biodegradable. Primarily used in oral patient care and examinations.

Dimensions

6" Overall Length (152.4 mm)
0.688" Width (17.475 mm)

Qualities

- Medical grade quality
- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges
- Made in USA

Applications

- Oral patient care and examinations
- Veterinary examinations
- Applying or mixing medicines
- Health and beauty applications - waxing
- Stirring stick
- Applying epoxy or resins

Ref #	Packaging
709 Non-sterile	10/250 (10 boxes per case, 250 per box)



Puritan®