

710

Junior Tongue Depressor

Sterile, junior size wood blade, individually wrapped

DESCRIPTION

This junior size tongue blade is made of Northern Maine white birch which is naturally biodegradable. Individually wrapped with no flavoring. Also available in three flavors: cherry (710 CHERRY), strawberry (710 STRAWBERRY), and grape (710 GRAPE). Primarily used in oral patient care and as a first aid kit component.

FEATURES & BENEFITS

- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges
- Puritan brand
- Made in USA

APPLICATIONS

- Oral patient care and examinations
- First aid kit component
- Health and beauty; waxing
- Stirring stick
- Applying epoxy or resins

Item No.	Dimensions length x diameter	Packaging pack, weight
non-sterile sterile		
25-711	139.7 x 15.875 mm 5.5 x 0.625"	10/250,* 21 lbs

*10 boxes per case, 250 per box (individually wrapped)

Products may not appear at actual size. See chart or our website for accurate measurements.

