



Foam Tipped Applicator - Dry Transport System

Standard foam tip, polystyrene handle

25-1506 1 PF BT BL

25-1506 1PF BT



Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kits. This specialty tipped applicator is made with 100ppi medical grade polyurethane foam. It is bonded securely to a tapered acetal handle.

Dimensions

Overall Length

- Tube: 6.062" (154 mm)
- Swab: 6.00" (152.4 mm)

Tip

- 0.188" Dia. (4.775 mm)
- 0.63" Length (16 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Qualities

- Medical grade quality
- Foam has excellent collection and release properties
- Soft, small tip for patient comfort
- Lint and residue free
- Made in USA

Applications

- Specimen collection
- Diagnostic sampling
- Recovery of microorganisms and fluids

Ref #	Packaging
25-1506 1PF BT Sterile	5/100/1 (5 boxes per case, 100 tubes per box, 1 swab per tube)
25-1506 1 PF BT BL Sterile	5/100/1 (5 boxes per case, 100 tubes per box, 1 swab per tube)



Puritan®