



# Tongue Depressor - Economy Grade

Standard size



- Biodegradable wood handle

## Description

This economy grade, standard size tongue depressor is made from white birch. This material has a high tensile strength and is naturally biodegradable with splinter free edges. This tongue blade is non-sterile and is primarily used in oral patient care and veterinary examinations.

## Dimensions

6" Overall Length (152.4 mm)

0.688" Width (17.46 mm)

## Qualities

- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges

## Applications

- Veterinary examinations
- Oral patient care

Ref #	Packaging
704 HOSPITAL Non-sterile	10/500 (10 boxes per case, 500 per box)



**HOSPITAL®**