

710 STRAWBERRY

Flavored Junior Tongue Depressor

Sterile, junior size wood blade with strawberry flavoring, individually wrapped

DESCRIPTION

This junior size tongue blade is made of Northern Maine white birch which is naturally biodegradable. Individually wrapped and strawberry flavored. Also available in plain (710), grape (710 GRAPE), and cherry (710 CHERRY) flavors. Primarily used for pediatrics, for oral patient care and examination. It can also be used as a first aid kit component.

FEATURES & BENEFITS

- Kosher certified, FDA food grade, allergen-free flavoring
- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges
- Individually wrapped
- Made in USA

APPLICATIONS

- Oral patient care and examinations
- First aid kit component
- Pediatrics
- Speech and occupational therapy
- Veterinary examinations

Item No.	Dimensions length x diameter	Packaging pack, weight
non-sterile sterile		
710 STRAWBERRY	139.7 x 15.875 mm 5.5 x 0.625"	20/50*, 10 lbs

*20 boxes per case, 50 per box (individually wrapped)

Products may not appear at actual size. See chart or our website for accurate measurements.

