




Flavored Junior Tongue Depressor

Wood



Description

This junior size tongue blade is made of Northern Maine white birch which is naturally biodegradable. Primarily used for pediatrics, for oral patient care and as a component in first aid kits. Individually wrapped and flavored (cherry, grape or strawberry).

Dimensions

5.5" Overall Length (139.7 mm)
0.625" Width (15.875 mm)

Qualities

- Medical grade quality
- Kosher certified, FDA food grade, Allergen free flavoring
- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges
- Made in USA

Applications

- Oral patient care and examinations
- First aid kit component
- Pediatrics
- Speech and occupational therapy
- Applying or mixing medicines
- Veterinary examinations

Ref #	Packaging
710 STRAWBERRY Non-sterile	20/50 (20 boxes per case, 50 per box)



Puritan®