



# Applicator Stick

## Wood



### Description

This wood applicator is made of Northern Maine White birch which is naturally biodegradable with a splinter free finish.

### Dimensions

12" Overall Length (304.8 mm)

0.125" Dia. (3.175 mm)

### Qualities

- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges
- Made in USA

### Applications

- Multi-use applicator stick
- For applying or stirring medicines
- Can be used in any setting: Medical, Health & Beauty, Industrial, Tattoo, Veterinary

| Ref #              | Packaging                                  |
|--------------------|--------------------------------------------|
| 809<br>Non-Sterile | 12/500<br>(12 boxes per case, 500 per box) |
|                    |                                            |



**Puritan®**