



# Applicator Stick

Wood



## Description

This wood applicator is made of Northern Maine White birch which is naturally biodegradable with straight cut ends. Made for general purpose functionality.

## Dimensions

- 18" Overall Length (457.2 mm)
- 0.125" Dia. (3.175 mm)

## Qualities

- Biodegradable handle
- Straight cut ends
- Birch has high tensile strength
- Made in USA

## Applications

- Multi-use applicator stick
- Applying or stirring medicines
- Specimen transfer
- Can be used in any setting: Medical, Health & Beauty, Industrial, Veterinary, Tattoo

Ref #	Packaging
811 Non-Sterile	2/1000 (2 bundles per case, 1000 per bundle)



**Puritan**<sup>®</sup>