



# Applicator Stick - Economy Grade

Round wood handle



- Biodegradable wood handle

## Description

This economy Applicator Stick is a sterile applicator stick that is made of Northern Maine white birch wood. This high-tensile strength material is naturally biodegradable and has a splinter-free finish making it perfect for a wide range of applications and environments.

## Dimensions

5.812 Overall Length (147.625 mm)  
0.083" Width (2.108 mm)

## Qualities

- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges

## Applications

- Patient and wound care
- First aid kit component
- Health and beauty - waxing or applying salves
- Removing particulates and residues from machinery or surfaces
- General cleaning
- Veterinary examinations

Ref #	Packaging
807 HOSPITAL Non-sterile	20/1,000 (20 boxes, 1,000 per box)



**HOSPITAL®**