




Cuticle/Orange Stick

Wood, double bevel both ends, wrapped



*Image not actual size

Description

This contaminate free cuticle stick is manufactured in a controlled environment to minimize contamination. Rigid white birch handle is double beveled on both ends. Individually wrapped in a paper/poly pouch.

Dimensions

7" x 0.156" (177.8 x 3.962 mm)

Qualities

- Double bevel on both sides
- Wood is naturally static dissipative
- Individually wrapped to minimize contact during production
- Manufactured in a controlled environment
- Made in USA

Applications

- Medical/Forensics/Health&Beauty
 - Evidence collection
 - Manicures and pedicures
 - Podiatry practices

Ref #	Packaging
2911 Non-sterile - wrapped	4/125/1 (4 boxes per case, 125 packages per box, 1 stick per package)



Puritan®