



# Applicator Stick

Tapered pointed applicator stick, sterile or non-sterile



\*Image not to scale

## Description

Pointed applicator stick for use as a colony transfer device with the automated MALDI-TOF MS systems. This stick is made out of Northern Maine White Birch. Available sterile and non-sterile.

## Dimensions

- Length: 2.875" (73 mm)
- Width: 0.083" (2.108 mm)

## Qualities

- Made out of Northern Maine White Birch
- Tapered pointed tip
- Made in USA

## Applications

- For use with MALDI-TOF MS
- Used as a direct bacterial colony transfer device to the MALDI racks for analysis



**Puritan**<sup>®</sup>

Ref #	Packaging
25-28107 Sterile	10/40/26 (10 boxes per case, 40 wrappers per box, 26 sticks per wrapper)
28107 PTD 5/600 Non-Sterile	5/600 (5 boxes per case, 600 sticks per box)