




PurFlock Ultra® Flocked Swab

Elongated Tip with Polystyrene Handle - 80 mm breakpoint



Description

PurFlock's purity make it ideal for PCR, molecular assays, DNA testing, forensic applications as well as direct antigen testing. Its uniquely designed fiber structure accounts for better cell yield in the specimen where sample size is limited.

Patented: US, Canada, Europe, Australia, Japan, China
Patent pending: Hong Kong

Dimensions

5.93" Overall Length (150.622 mm)

Tip

- 0.219" Dia. (5.563 mm)
- 0.770" Length (19.558 mm)

Handle

- 0.1" Dia. (2.54 mm) tapered to 0.09" (2.286 mm)
- 5.89" Length (149.606 mm)

Qualities

- Medical grade quality
- Highly absorbent tip
- Soft smooth tip
- Made in USA

Applications

- For Buccal cell collection
- Specimen collection
- Rapid diagnostic testing
- Forensic evidence collection
- DNA testing



Ref #	Packaging
3506-U Non-Sterile	10/100 (10 bags per case, 100 per bag)
25-3506-U Sterile	10/50/1 (10 boxes of 50 individually wrapped swabs)