



Cotton Tipped Applicator - Dry Transport System

Standard cotton tip, polystyrene or wooden handle



Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kits. This absorbent cotton tipped applicator is produced with high-quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a semi-flexible polystyrene or a rigid wood handle by an aqueous based adhesive.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.00" (152.4 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Qualities

- High-quality USP grade cotton
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

Applications

- Specimen collection for testing and screening
- Cell collection for DNA testing

Ref #	Packaging
25-806 1PC BT Sterile	5/100/1 (5 boxes per case, 100 tubes per box, 1 swab per tube)
25-806 1WC BT Sterile	5/100/1 (5 boxes per case, 100 tubes per box, 1 swab per tube)



Puritan®