



Foam Tipped Applicator- Dry Transport System

Standard tip, polystyrene handle, molded cap



Molded base

Description

Puritan's traditional dry transport tube is perfect for specimen collection and transportation to the lab or use in diagnostic kits. This specialty applicator is tipped with 100ppi medical grade polyurethane foam. It is bonded securely to a polystyrene handle. The molded base holds the swab snugly in the center of the tube.

Overall Length

- Tube: 6.625" (168.275 mm)
- Swab: 5.979" (151.86 mm)

Tip

- 0.188" Dia. (4.77 mm)
- 0.655" Length (16.63 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Qualities

- Absorbent tip
- High particle collection capacity
- Soft and smooth tip
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Molded base holds swab in center of the tube
- Made in USA

Applications

- Ideal for specimen collection and rapid diagnostics and transport forensics
- DNA testing of family / species / gender
- Diagnostic specimen collection
- Wildlife forensic biology
- Buccal and other cell collection for testing and screening



Puritan®

Ref #	Packaging
25-1506 1PF TT MC Sterile	10/100/1 (10 boxes of 100 tubes per box, 1 swab per tube = 1,000/cs)