



Cotton Tipped Applicator

Extra Absorbent Large Tip, Wood Handle, Bulk Packaging

Description

This extra-absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a wood handle by an aqueous based adhesive. Available with 10 or 50 swabs per package.

Dimensions

6" Overall Length (76.2 mm)

Tip

- 0.234" Dia. (5.94 mm)
- 0.687" Length (17.449 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.812" Length (147.624 mm)

10 swabs per package



Qualities

- Medical grade quality
- Extra absorbent, large tip
- Soft smooth tip
- Rigid wood handle
- 10 or 50 swabs per package
- Made in USA

Applications

Medical/Critical Environments & Industrial

- Applying salves or ointments
- Wound care/dressing and cleaning
- Specimen collection
- Site preparation
- Cosmetic application/tattoos
- Medication applicator



50 swabs per package

Ref #	Packaging
25-806 10WCL Sterile	10/5/50 (10 boxes per case, 5 bags per box, 50 swabs per bag)
25-806 WCLB 10/50 Sterile	10/50 (10 bags of 50 swabs per bag = 500/cs)



Puritan®