



HydraFlock® Flocked Swab

Straight Tip with Extruded Polystyrene Handle - 30mm breakpoint



Description

A new era in flock swab technology. HydraFlock provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Patented: US, Canada, Europe, Australia, Japan, China
Patent pending: Hong Kong

Dimensions

5.92" Overall Length (150.36 mm)

Tip

- 0.18" Dia. (4.572 mm)
- 0.800" Length (20.32 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.88" Length (149.35 mm)

Qualities

- Highly absorbent tip
- Soft smooth tip
- Straight tip
- Extruded polystyrene handle with 30mm breakpoint
- Made in USA

Applications

- Buccal cell collection/DNA testing
- Specimen collection
- Forensic evidence collection
- Rapid diagnostics testing

Ref #	Packaging
3000-H E30 Non-Sterile	10/100 (10 bags per case, 100 swabs per bag = 1,000/cs)
25-3000-H E30 Sterile	10/50/1 (10 boxes per case of 50 individually wrapped swabs = 500/cs)



HydraFlock®