




PurFlock Ultra® Flocked Swab

Straight Tip with Extruded Polystyrene Handle - 30mm breakpoint



Description

An ideal choice in flock swabs with excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Patented: US, Canada, Europe, Australia, Japan, China
Patent pending: Hong Kong

Dimensions

5.92" Overall Length (150.36 mm)

Tip

- 0.18" Dia. (4.572 mm)
- 0.800" Length (20.32 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.88" Length (149.35 mm)

Qualities

- Highly absorbent tip
- Soft smooth tip
- Straight tip
- Extruded polystyrene handle with 30mm breakpoint
- Made in USA

Applications

- Buccal cell collection/DNA testing
- Specimen collection
- Forensic evidence collection
- Rapid diagnostics testing

Ref #	Packaging
3000-U E30 Non-Sterile	10/100 (10 bags per case, 100 swabs per bag = 1,000/cs)
25-3000-U E30 Sterile	10/50/1 (10 boxes per case of 50 individually wrapped swabs, 1 swab p/pkg = 500/cs)



**PurFlock
ULTRA®**