



Cotton Tipped Applicator

Extra large tip, wood handle



Description

This extra-absorbent, extra large cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a sturdy wood handle with a hot melt glue adhesive.

Dimensions

5.937" Overall Length (150.799 mm)

Tip

- 0.312" Dia. (7.924 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.812" Length (147.62 mm)

Qualities

- Extra large cotton tip
- Naturally biodegradable wood handle
- Ridig handle
- Made in USA

Applications

- General cleaning swab
- Applying medicines
- Cleaning medicine dispensers
- Removing excess from injection molding
- Removing dirt and oil from machine parts

Ref #	Packaging
806-WCXL Non-sterile	10/10/50 (10 boxes per case, 10 bags per box, 50 swabs per bag)



Puritan®