

Microbiology Media Products

Cross Reference Sheet



RELATED PRODUCTS							PRODUCT DESCRIPTION			
Puritan	Copan	BD	Hardy Diagnostics	Medical Wire	Remel: ThermoFisher Sci	Zymo	Description	Vial Diameter x Length (mm)	Fill Volume	Cap Color
CB-200	n/a	n/a	n/a	n/a	n/a	n/a	Cary-Blair: Transport medium for the preservation of fecal and rectal specimens containing enteric pathogenic bacteria. For <i>in vitro</i> diagnostic use.	12 x 80	2mL	orange
CB-500	n/a	n/a	n/a	n/a	R060450/ R060452 (9mL)	n/a		16 x 100	5mL	
CT-200	094CE.A	n/a	n/a	n/a	n/a	n/a	CT Broth: Selective enrichment broth medium for use in the isolation of <i>Candida spp</i> and <i>Trichomonas spp.</i> For <i>in vitro</i> diagnostic use.	12 x 80	2mL	blue
CT-500	092CN (4ml)	298323	n/a	n/a	n/a	n/a		16 x 100	5mL	
LB-200	476CE.A	n/a	n/a	MWGBST	R064812 (3mL)	n/a	Lim Broth: Selective enrichment medium for use in the qualitative procedures for the isolation of Group B Streptococcus (GBS) from clinical specimens. For <i>in vitro</i> diagnostic use.	12 x 80	2mL	light blue
LB-500	n/a	292209/ 296266	L57	n/a	R064814	n/a		16 x 100	5mL	
MRSA-200	477CE.A (TSB Salt Broth)	n/a	n/a	n/a	n/a	n/a	MRSA: Enrichment broth medium used for the isolation of specimens containing MRSA. For <i>in vitro</i> diagnostic use.	12 x 80	2mL	pink
MRSA-500	n/a	n/a	n/a	n/a	n/a	n/a		16 x 100	5mL	
RL-200	n/a	n/a	n/a	n/a	n/a	n/a	Regan-Lowe: Semi-solid transport medium for the isolation of <i>Bordetella</i> species. For <i>in vitro</i> diagnostic use.	12 x 80	2mL	white
RL-500	n/a	297855 (10mL)	Q32 (4mL; 13 x 100mm tube)	n/a	R064141	n/a		16 x 100	5mL	
SB-200	475CE.A	n/a	n/a	MWSEL	n/a	n/a	Selenite Broth: Enrichment broth medium used for the isolation of <i>Salmonella spp.</i> and <i>Shigella spp.</i> For <i>in vitro</i> diagnostic use.	12 x 80	2mL	white
SB-500	n/a	221020 (8mL)	n/a	n/a	R064500 / R064502 (9mL)	n/a		16 x 100	5mL	

MADE IN THE U.S.A

6.17 VAT



twitter.com/PuritanUSA



facebook.com/PuritanUSA



linkedin.com/company/PuritanUSA



youtube.com/PuritanUSA

800-321-2313 • puritanmedproducts.com

SP-100	097CE.A	n/a	n/a	n/a	n/a	n/a	Sputum Solution: Liquefying solution of clinical sputum specimens for the isolation of bacteria and fungi that cause respiratory tract infections. For <i>in vitro</i> diagnostic use.	12 x 80	1mL	natural
SP-100P	098CE.A	n/a	n/a	n/a	n/a	n/a		12 x 80	1mL	
VL 1201 SAFE	608C	n/a	n/a	n/a	n/a	R1102	PurSafe® DNA/RNA Preservative: Transport medium for the preservation of DNA and RNA of clinical specimens. Not for diagnostic use. Research use only.	12 x 80	1mL	light blue
VL 1201 SAFE-H	608CS01R	n/a	n/a	n/a	n/a	R1106 / R1107		12 x 80	1mL	
VL 1202 SAFE	606C	n/a	n/a	n/a	n/a	n/a		12 x 80	2mL	light blue
VL 1202 SAFE-H	606CS01R	n/a	n/a	n/a	n/a	R1108/ R1109		12 x 80	2mL	