

PURITAN MEDICAL PRODUCTS COMPANY, LLC

Safety Data Sheet

November 16, 2020

1. Product and Company Identification

Product Name: Fecal Opti-Swab Collection and Transport System

Product Number: CB-206

Company Identification:

Puritan Medical Products Company LLC
P.O. Box 149, 31 School Street
Guilford, Maine 04443-0149 U.S.A

Contact numbers:

Tel: +1 207-876-3311
Fax: +1 207-876-3130

2. Hazardous identification

Skin contact: Irritation

Eye contact: Irritation

3. Composition/information on ingredients

Description: Liquid Mixture

Hazardous Components:	Concentration:	R	Classification:
Sodium Chloride	0.3%	36/38	7647-14-5
Calcium Chloride	0.01%	36	10043-52-4
Sodium Thioglycolate	0.15%	22/36/37/38	367-51-1

4. First-aid measures

Skin contact: Immediately wash with water and soap, then rinse thoroughly.

Eye contact: Flush open eye with plenty of water for 15 minutes. Consult doctor.

Inhalation: Supply fresh air; if breathing is difficult give oxygen. Get medical attention if cough or other symptoms appear.

Swallowing: Wash mouth with water, if symptoms persist consult doctor.

5. Fire-fighting measures

Suitable Extinguishing media: Dry chemicals, CO₂, foam or Halon extinguisher

Protective Equipment: Personal Protective Equipment (PPE)

6. Accidental release measures

Personal precautions: Personal Protective Equipment (PPE)

Environmental precautions: N/A

Methods for Clean-up: Wipe up with damp cloth or mop.

7. Handling and storage

Handling: No special handling procedures required

Storage: Store away from oxidizing agents. Store at 2°- 25°C.

8. Exposure controls/personal protection

Respiratory protection: N/A

Hand protection: Disposable rubber gloves

Eye protection: Safety glasses

Skin and body protection: Protective work clothing (lab coat)

Safety Data Sheet

9. Physical and chemical properties	
Appearance/Odor	Characteristic
pH:	7.2 ± 0.3 @ 25°C
Density:	Not determined
Boiling point, °C:	Not determined
Melting point, °C	Not determined
Flash point, °C:	Not applicable
Solubility/Miscibility with water:	N/A
10. Stability and reactivity	
Thermal decomposition/conditions to avoid:	No decomposition and mucous membranes
Incompatible materials:	N/A
Hazardous decomposition products:	No dangerous decomposition products known
11. Toxicology information	
Skin	Irritant to skin and mucous membranes
Eye	Irritant effect
Sensitization	No sensitization effect known
Additional Toxicological Information	The product shows the following dangers according to internal approved preparations: irritant
12. Ecological Information	
Ecology Effects:	None have been identified.
General Notes:	Not known to be hazardous to water.
13. Disposal considerations	
Recommendation:	Dispose of materials in accordance with federal (40 CFR 261.3), state, and local requirements. This product is not considered a RCRA hazardous waste.
14. Transport information	
ADR/RID:	No
IMO/IMDG:	No
KAO/IATA:	No
15. Regulatory information	
SARA Section 355:	None of the ingredients is listed
SARA Section 313:	None of the ingredients is listed
TSCA:	7647-14-5 Sodium Chloride 10043-52-4 Calcium Chloride 367-51-1 Sodium Thioglycolate
California Proposition 65:	None of the ingredients is listed
16. Other information	
Puritan Medical Products Company LLC provides the information in this document in good faith and believes the information to be accurate. The chemical, physical and toxicological properties of this product have not been thoroughly investigated. It is the responsibility of the buyer to research and understand safe methods of handling, storing, and disposal of this product. Puritan Medical makes no warranty with respect to such information and assumes no liability for any loss or injury, which may result from the use of this information. It is the buyer's responsibility to comply with local, state, and federal regulations concerning use and disposal of this product.	