



Foam Tipped Applicator - Dry Transport System

Double Standard Foam Tip, Polystyrene Handle



Description

Puritan's polypropylene double dry transport tube and plug-style cap is perfect for specimen collection and transport to the lab or for use in diagnostic kits. These specialty tipped applicators are made with 100ppi medical grade polyurethane foam. They are bonded securely to polystyrene handles.

Dimensions

Overall Length

- Tube: 6.062" (154 mm)
- Swab: 5.979" (151.9 mm)

Tip

- 0.188" Dia. (4.8 mm)
- 0.630" Length (16 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.875" Length (149.2 mm)

Qualities

- Medical grade quality
- Foam has excellent collection and release properties
- Double foam tip
- Soft, small tips for patient comfort
- Lint and residue free
- Made in USA

Applications

- Specimen collection
- Diagnostic sampling
- Recovery of microorganisms and fluids

Ref #	Packaging
25-1506 2PF BT Sterile	10/50/2 (10 boxes per case, 50 tubes per box, 2 swabs per tube)



Puritan®