

# PURITAN MEDICAL PRODUCTS COMPANY, LLC

## Safety Data Sheet

30-Dec-16

### 1. Product and Company Identification

**Product Name:** Sputum Solution

**Product Number:** SP-100P

**Company Identification:**

Puritan Medical Products Company LLC  
P.O. Box 149, 31 School Street  
Guilford, Maine 04443-0149 U.S.A

**Contact numbers:**

Tel: +1 207-876-3311  
Fax: +1 207-876-3130

### 2. Composition/information on ingredients

Description: Liquid Mixture

Hazardous Components: No ingredient is hazardous

### 3. Hazards Identification

Ingestion: N/A

Skin contact: N/A

### 4. First-aid measures

Skin contact: Immediately wash with water and soap, then rinse thoroughly.

Eye contact: Flush open eye with plenty of water for 15 minutes. Consult doctor.

Inhalation: Supply fresh air; if breathing is difficult give oxygen. Get medical attention if cough or other symptoms appear.

Swallowing: Immediately call a doctor.

### 5. Fire-fighting measures

Suitable Extinguishing media: Dry chemicals, CO2, foam or Halon extinguisher

Protective Equipment: Personal Protective Equipment (PPE)

### 6. Accidental release measures

Personal precautions: Personal Protective Equipment (PPE)

Environmental precautions: N/A

Methods for Clean-up: Wipe up with damp cloth or mop. Ensure proper ventilation.

### 7. Handling and storage

Handling: No special handling procedures required

Storage: Store away from oxidizing agents. Store at 2°- 25°C.

### 8. Exposure controls/personal protection

Respiratory protection: N/A

Hand protection: Disposable rubber gloves

Eye protection: Safety glasses

Skin and body protection: Protective work clothing (lab coat)

## Safety Data Sheet

### 9. Physical and chemical properties

Appearance/Odor	Characteristic
pH:	7.0 ± 0.5 @ 25°C
Density:	Not determined
Boiling point, °C:	Not determined
Melting point, °C	Not determined
Flash point, °C:	Not applicable
Solubility/Miscibility with water:	Soluble

### 10. Stability and reactivity

Thermal decomposition/conditions to avoid:	Stable under normal conditions.
Incompatible materials:	Strong oxidizing agents, strong bases
Hazardous decomposition products:	No data available

### 11. Toxicology information

Skin	No data available
Eye	No data available
Ingestion	No data available
Sensitization	No data available
Additional Toxicological Information	None

### 12. Ecological Information

Ecology Effects:	None have been identified.
General Notes:	None have been identified.

### 13. Disposal considerations

Recommendation:	Dispose of materials in accordance with federal (40 CFR 261.3), state, and local requirements. This product is not considered a RCRA hazardous waste.
-----------------	--

### 14. Transport information

ADR/RID:	No
IMO/IMDG:	No
KAO/IATA:	No

### 15. Regulatory information

SARA Section 355:	None of the ingredients are listed
SARA Section 313:	None of the ingredients are listed
TSCA:	None of the ingredients are listed
California Proposition 65:	None of the ingredients are listed

### 16. Other information

Puritan Medical Products Company LLC provides the information in this document in good faith and believes the information to be accurate. The chemical, physical and toxicological properties of this product have not been thoroughly investigated. It is the responsibility of the buyer to research and understand safe methods of handling, storing, and disposal of this product. Puritan Medical makes no warranty with respect to such information and assumes no liability for any loss or injury, which may result from the use of this information. It is the buyer's responsibility to comply with local, state, and federal regulations concerning use and disposal of this product.