



Foam Tipped Applicator

Small foam tip, polystyrene handle



Description

This foam tipped applicator is made with reticulated 100ppi medical grade foam. Foam has excellent collection and release properties.

Dimensions

6.0" Overall Length (152.4 mm)

Tip

- 0.078" Dia. (1.981 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.085" Dia. tapered to 0.040" (2.159 to 1.016 mm)
- 5.848" Length (148.539 mm)

Qualities

- Soft, non-abrasive tip for patient comfort
- Can be used with IPA or DI water, or dry
- Made in USA

Applications

- Buccal swab collection
- For specimen collection
- For use in a wide range of diagnostic kits
- For general cleaning in industrial applications

Ref #	Packaging
1506-PFB Non-sterile	4/250 (4 bags per case, 250 swabs per bag)



Puritan®