



# Calgiswab® Calcium Alginate Tipped Applicator

Nasopharyngeal, aluminum handle



## Description

This specialty applicator is tipped with calcium alginate fiber, an organic material derived from seaweed and carded in Puritan's facility to assure conformance with medical device specifications. Best known as a nasopharyngeal (NP) swab.

## Dimensions

**5.437" Overall Length (138.1 mm)**

### Tip

- 0.109" Dia. (2.77 mm)
- 0.437" Length (11.1 mm)

### Handle

- 0.0351" Dia. (0.89 mm)
- 5.375" Length (136.525 mm)

## Qualities

- Medical grade quality
- Calcium alginate tip contains no fatty acids, which may be inhibitory to many strains of bacteria
- Thin flexible aluminum handle for minimal patient discomfort
- Made in USA

## Applications

- Flexible nasopharyngeal(NP) swab
- Specimen collection for testing and screening
- Suitable for small cavities and wound care
- Pediatrics
- Ophthalmology
- Veterinary care

| Ref #                  | Packaging  |
|------------------------|--|
| 25-800 A 50<br>Sterile | 10/50/1<br>(10 boxes of 50 individually wrapped packages, 1 p/pkg) |



**Puritan®**