
Cotton Tipped Applicator

Wood handle





Description

This absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a rigid white birch handle by an aqueous based adhesive.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.188" Dia. (4.775 mm)
- 0.625" Length (15.875 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.812" Length (147.624 mm)

Qualities

- Medical grade quality
- Standard size cotton tip
- Absorbent tip
- · Soft smooth tip
- · Made in USA

Applications

- Specimen collection
- · Wound care
- · Cell collection for DNA testing
- · Medication applicator
- · First-aid kit component
- Site preparation for health/beauty and cosmetic applications
- Tattoos

| Ref # | Packaging |
|------------------------|--|
| 25-806 1 WC Sterile | 10/100/1 (10 boxes of 100 individually wrapped packages, 1 p/pkg) |
| 25-806 2WC Sterile | 10/100/2 (10 boxes of 100 individually wrapped packages, 2 p/pkg) |
| 25-806 10WC Sterile | 10/20/10 (10 boxes of 20 individually wrapped packages, 10 p/pkg) |



Tel 800 321-2313 (US & Canada Only) or 207-876-3311

Web puritanmedproducts.com Email sales@puritanmedproducts.com