



# PurFlock Ultra® Flocked Swab

Ultrafine Tip with Polystyrene Handle - 100mm Breakpoint



## Description

An ideal choice in flock swabs with excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Patented: US, Canada, Europe, Australia, Japan, China, Hong Kong

## Dimensions

**6.04" Overall Length (153.416 mm)**

### Tip

- 0.140" Dia. (3.556 mm)
- 0.66" Length (16.76 mm)

### Handle

- 0.1" Dia. (2.54 mm) tapered to 0.038" (0.965 mm)
- 6" Length (152.4 mm)

## Qualities

- Newly designed nasopharyngeal swab
- Medical grade quality
- Highly absorbent tip
- Soft smooth tip
- Made in USA

## Applications

- For nasopharyngeal sample collection
- Specimen collection
- Rapid diagnostic testing
- Forensic evidence collection

Ref #	Packaging
25-3320-U EMB 100MM Sterile	10/50/1 (10 boxes of 50 individually wrapped swabs)

