

PurFlock Ultra® Flocked Swab

Ultrafine Tip with Polystyrene Handle - 80mm Breakpoint



80mm breakpoint

Description

An ideal choice in flock swabs with excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Patented: US, Canada, Europe, Australia, Japan, China, Hong Kong

Dimensions

6.04" Overall Length (153.416 mm)

Tip

- 0.140" Dia. (3.556 mm)
- 0.66" Length (16.76 mm)

Handle

- 0.1" Dia. (2.54 mm) tapered to 0.038" (0.965 mm)
- 6" Length (152.4 mm)

Qualities

- · Newly designed nasopharyngeal swab
- Medical grade quality
- · Highly absorbent tip
- Soft smooth tip
- Made in USA

Applications

- For nasopharyngeal sample collection
- · Specimen collection
- · Rapid diagnostic testing
- · Forensic evidence collection

| Ref # | Packaging |
|-------------------------------|---|
| 25-3320-U EMB 80MM Sterile | 10/50/1 (10 boxes of 50 individually wrapped swabs) |



Tel 800 321-2313 (US & Canada Only) or 207-876-3311

Web puritanmedproducts.com Email sales@puritanmedproducts.com