




HydraFlock® Flocked Swab

Ultrafine Tip with Polystyrene Handle - 80mm breakpoint



Description

A new era in flock swab technology. HydraFlock provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Patented: US, Canada, Europe, Australia, Japan, China, Hong Kong

Dimensions

6.04" Overall Length (153.416 mm)

Tip

- 0.140" Dia. (3.556 mm)
- 0.66" Length (16.76 mm)

Handle

- 0.1" Dia. (2.54 mm) tapered to 0.038" (0.965 mm)
- 6" Length (152.4 mm)

Qualities

- Newly designed nasopharyngeal swab
- Medical grade quality
- Highly absorbent tip
- Soft smooth tip
- Made in USA

Applications

- For nasopharyngeal sample collection
- Specimen collection
- Rapid diagnostic testing
- Forensic evidence collection

Ref #	Packaging
25-3320-H EMB 80MM Sterile	10/50/1 (10 boxes of 50 individually wrapped swabs)



HydraFlock®