




PurFlock Ultra® Flocked Swab

Elongated Tip with Polystyrene Handle - 80 mm breakpoint



Description

PurFlock's purity make it ideal for PCR, molecular assays, DNA testing, forensic applications as well as direct antigen testing. Its uniquely designed fiber structure accounts for better cell yield in the specimen where sample size is limited.

Patented: US, Canada, Europe, Australia, Japan, China, Hong Kong

Dimensions

5.93" Overall Length (150.622 mm)

Tip

- 0.219" Dia. (5.563 mm)
- 0.770" Length (19.558 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.89" Length (149.606 mm)

Qualities

- Medical grade quality
- Highly absorbent tip
- Soft smooth tip
- Made in USA

Applications

- For Buccal cell collection
- Specimen collection
- Rapid diagnostic testing
- Forensic evidence collection
- DNA testing

Ref #	Packaging
25-3506-U Sterile	10/50/1 (10 boxes of 50 individually wrapped swabs)

