



# Norovirus Detection Swab

Round Macrofoam Tip, Polypropylene Handle, Transport Tube



## Description

Ready-to-use norovirus detection swab and collection tube. Available with a round macrofoam tip adhered to a polypropylene handle, and packaged in a dry transport tube. Choose from a pre-moistened or dry swab depending on your sampling site. Used for the detection of human norovirus on environmental surfaces.

## Dimensions

### Overall Length

- Tube: 6.02" (152.908 mm)
- Swab: 5.938" (150.82 mm)

### Tip

- 0.850" Dia. (21.59 mm)
- 1.375" Length (34.92 mm)

### Handle

- .128" Dia. (3.25 mm)
- 5.81" Length (147.57 mm)

## Qualities

- Large tip for maximized collection and elution
- Quick-turn-open leak-proof cap
- Swab securely attached to cap
- High-quality Puritan Foam tipped swab
- Available pre-moistened or dry
- Unwrapped
- Shelf-life: 6 months (pre-moistened); 3 years (dry)
- Made in USA

## Applications

- Used for the detection of human norovirus on environmental surfaces using a newly developed protocol by the Center for Disease Control
- Assists in determining the level of environmental contamination during outbreaks
- Detects viruses when clinical samples are not available and monitors the effectiveness of cleaning practices

Ref #	Packaging
25-88060 PF UW Sterile	25/1 (25 individually wrapped tubes, 1 swab per tube = 25/cs)
25-88060 PF UW DRY Sterile	



**Puritan®**