




Foam Tipped Applicator

Polystyrene handle



Description

This specialty applicator is tipped with 100ppi medical grade polyurethane foam. It is bonded securely to a hollow polystyrene handle that snaps when bent at a 90° angle.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.188" Dia. (4.775 mm)
- 0.63" Length (16.002 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Qualities

- Medical grade quality
- Foam has excellent collection and release properties
- Soft foam tip for patient comfort
- Can be used dry or with media
- Made in USA

Applications

- Buccal swab
- For nasal and anterior nares specimen collection for diagnostic testing
- For general specimen collection for testing and screening
- Used in Pediatrics, Ophthalmology, Veterinary care, and Oropharyngeal applications
- Component in rapid diagnostic kits

Ref #	Packaging
25-1506 1PF 100 Sterile	10/100/1 paper/paper pouch (10 boxes of 100 individually wrapped packages per case, 1 swab per/pkg)



Puritan®