




Puritan® Foam Tipped Applicator

Standard Foam Tip, Breakable Hollow Polystyrene handle



Description

This specialty applicator is tipped with 100ppi medical grade polyurethane foam. It is bonded securely to a hollow polystyrene handle that snaps when bent at a 90° angle. Available in multiple packaging configurations.

Dimensions

5.979" Overall Length (151.867 mm)

Tip

- 0.188" Dia. (4.775 mm)
- 0.630" Length (16.002 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Qualities

- Medical grade quality
- Foam has excellent collection and release properties
- Soft foam tip for patient comfort
- Can be used dry or with media
- Made in USA

Applications

- Buccal swab
- For nasal and anterior nares specimen collection for diagnostic testing
- For general specimen collection for testing and screening
- Used in Pediatrics, Ophthalmology, Veterinary care, and Oropharyngeal applications
- Component in rapid diagnostic kits

Ref #	Packaging
25-1506 1PF 100 Sterile	(10/100/1) paper/paper pouch (10 boxes of 100 individually wrapped swabs, 1 swab p/pkg)
25-1506 1PFB 20 Sterile	1,000 BULK (50/20/1) paper/poly pouch (50 perforated strips, 20 swabs per strip, 1 swab p/pkg)



Puritan®