



Cotton Tipped Applicator

Small tip, Wood handle



Description

This absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. Special construction eliminates chemical migration.

Dimensions

3.032" Overall Length (77.013 mm)

Tip

- 0.156" Dia. (4.775 mm)
- 0.625" Length (15.875 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 2.906" Length (73.812 mm)

Qualities

- Small cotton tip, no glue adhesive
- Made in USA

Applications

Critical Environments / Industrial / Medical

- Used for cleaning and removing residues
- Ideal as a disk drive head cleaner
- Can be used for cleaning wounds and applying medicines

Ref #	Packaging
838-WCS NO GLUE Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)



Puritan®