

700

# Infant Tongue Depressor

Non-sterile, infant size wood blade

## DESCRIPTION

This infant size blade is made of Northern Maine white birch which is naturally biodegradable. Also available sterile (25-700). Primarily used in oral patient care on infants and small children.

## FEATURES & BENEFITS

- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges
- Made in USA

## APPLICATIONS

- Oral patient care and examinations
- Veterinary examinations
- Health and beauty; waxing
- Stirring stick

Item No.	Dimensions length x diameter	Packaging pack, weight
non-sterile   sterile		
700	114.3 x 9.525 mm 4.5 x 0.375"	10/1000,* 37 lbs

\*10 boxes per case, 1,000 per box

Products may not appear at actual size. See chart or our website for accurate measurements.

