




# Polyester Tipped Applicator

Plastic handle, inherent score mark



## Description

This synthetic tipped applicator is made of medical-grade polyester fiber that is carded in our factory to maintain consistency. It is non-absorbent and contains no fatty acids or latent sugars, making it ideal for specimen transport, lab work, and diagnostic kits. It is bonded securely to a tapered plastic handle by an aqueous based adhesive. Inherent score mark on handle. Available sterile or non-sterile.

## Dimensions

**5.867" Overall Length (149.022 mm)**

### Tip

- 0.125" Dia. (3.175 mm)
- 0.687" Length (17.45 mm)

### Handle

- 0.097" Dia. (2.4638 mm)
- 5.805" Length (147.447mm)

## Qualities

- Medical grade quality
- Non-absorbent tip
- Synthetic polyester contains no fatty acids or latent sugars
- Small tip for patient comfort
- Molded handle is semi-flexible
- Made in USA

## Applications

- Diagnostic specimen collection
- Recovery of micro-organisms and fluid
- Cell collection for DNA testing
- For use in a wide range of diagnostic tests and laboratory uses

Ref #	Packaging
1000-PD Non-Sterile	10/100 (10 bags per case, 100 swabs per bag)
25-1000 1PD Sterile	10/100/1 (10 boxes of 100 individually wrapped packages, 1 swab p/pkg)



**Puritan®**