

## 704 / 704 BULK

# Puritan® Standard Tongue Depressor

Non-sterile, standard size wood blade, available in bulk

### DESCRIPTION

This standard size blade is made of North American hardwood which is naturally biodegradable. Used in oral patient care and examinations, health and beauty applications, and a variety of industrial uses.

### FEATURES & BENEFITS

- Biodegradable handle
- High tensile strength
- Splinter free edges
- #704 – Puritan brand
- Made in USA

### APPLICATIONS

- Oral patient care and examinations
- Veterinary examinations
- Health and beauty; waxing
- Stirring stick
- Applying epoxy or resins
- Industrial use

Item No.	Dimensions length x diameter	Packaging pack, weight
non-sterile   <b>sterile</b>		
704	152.4 x 17.475 mm 6 x 0.688"	10/500,* 35 lbs
704 BULK		3,700 38 lbs (bulk)

\*10 boxes per case, 500 per box

Products may not appear at actual size. See chart or our website for accurate measurements.

