

# 705

## Puritan® Adult Tongue Depressor

Non-sterile, adult size wood blade

### DESCRIPTION

This adult size tongue blade is made of North American hardwood which is naturally biodegradable. Also available sterile (25-705). Primarily used in oral patient care and examinations.

### FEATURES & BENEFITS

- Biodegradable handle
- High tensile strength
- Splinter free edges
- Made in USA

### APPLICATIONS

- Oral patient care and examinations
- Veterinary examinations
- Health and beauty; waxing
- Stirring stick
- Applying epoxy or resins

Item No.	Dimensions length x diameter	Packaging pack, weight
non-sterile   sterile		
705	152.4 x 19.05 mm 6 x 0.75"	10/500,* 38 lbs

\*10 boxes per case, 500 per box

Products may not appear at actual size. See chart or our website for accurate measurements.

