

710 CHERRY

Puritan® Flavored Junior Tongue Depressor

Sterile, junior size wood blade with cherry flavoring, individually wrapped

DESCRIPTION

This junior size tongue blade is made of North American hardwood which is naturally biodegradable. Individually wrapped and cherry flavored. Also available in plain (710), strawberry (710 STRAWBERRY), and grape (710 GRAPE) flavors. Primarily used for pediatrics, for oral patient care and examination. It can also be used as a first aid kit component.

FEATURES & BENEFITS

- Kosher certified, FDA food grade, allergen-free flavoring
- Biodegradable handle
- High tensile strength
- Splinter free edges
- Individually wrapped
- Made in USA

APPLICATIONS

- Oral patient care and examinations
- First aid kit component
- Pediatrics
- Speech and occupational therapy
- Veterinary examinations

| Item No. | Dimensions length x diameter | Packaging pack, weight |
|-----------------------|-----------------------------------|---------------------------|
| non-sterile sterile | | |
| 710 CHERRY | 139.7 x 15.875 mm 5.5 x 0.625" | 20/50,* 10 lbs |

*20 boxes per case, 50 per box (individually wrapped)

Products may not appear at actual size. See chart or our website for accurate measurements.

