

## 710 STRAWBERRY

# Puritan® Flavored Junior Tongue Depressor

Sterile, junior size wood blade with strawberry flavoring, individually wrapped

### DESCRIPTION

This junior size tongue blade is made of North American hardwood which is naturally biodegradable. Individually wrapped and strawberry flavored. Also available in plain (710), grape (710 GRAPE), and cherry (710 CHERRY) flavors. Primarily used for pediatrics, for oral patient care and examination. It can also be used as a first aid kit component.

### FEATURES & BENEFITS

- Kosher certified, FDA food grade, allergen-free flavoring
- Biodegradable handle
- High tensile strength
- Splinter free edges
- Individually wrapped
- Made in USA

### APPLICATIONS

- Oral patient care and examinations
- First aid kit component
- Pediatrics
- Speech and occupational therapy
- Veterinary examinations

Item No.	Dimensions length x diameter	Packaging pack, weight
non-sterile   sterile		
710 STRAWBERRY	139.7 x 15.875 mm 5.5 x 0.625"	20/50*, 10 lbs

\*20 boxes per case, 50 per box (individually wrapped)

Products may not appear at actual size. See chart or our website for accurate measurements.

