

712

Puritan® Junior Tongue Depressor

Non-sterile, junior size wood blade

DESCRIPTION

This junior size tongue blade is made of North American hardwood which is naturally biodegradable. Primarily used in oral patient care and as a first aid kit component.

FEATURES & BENEFITS

- Biodegradable handle
- High tensile strength
- Splinter free edges
- Made in USA

APPLICATIONS

- Oral patient care and examinations
- First aid kit components
- Veterinary examinations
- Health and beauty; waxing

Item No.	Dimensions length x diameter	Packaging pack, weight
non-sterile sterile		
712	139.7 x 15.875 mm 5.5 x 0.625"	12/100,* 8 lbs

*12 boxes per case, 100 per box

Products may not appear at actual size. See chart or our website for accurate measurements.



CALL 800-321-2313 (USA & Canada) • 207-876-3311

EMAIL sales@puritanmedproducts.com

VISIT puritanmedproducts.com

PURITAN MEDICAL PRODUCTS CO., GUILFORD, ME, 04443, USA ISO 13485

712 04/13/22