

807

## Puritan® Applicator Stick

Non-sterile, wood stick

## DESCRIPTION

This wood applicator is made of North American hardwood which is naturally biodegradable with a splinter free finish. Also available in sterile two-pack (25-807 2W).

## FEATURES &amp; BENEFITS

- Biodegradable handle
- High tensile strength
- Splinter free edges
- Made in USA

## APPLICATIONS

- Multi-use applicator stick
- Can be used in any setting: medical, health & beauty, industrial, tattoo, veterinary

Item No.	Dimensions length x diameter	Packaging pack, weight
non-sterile   sterile		
807	147.625 x 2.108 mm 5.812 x 0.083"	20/1000* 18 lbs

\*20 boxes per case, 1,000 per box

Products may not appear at actual size. See chart or our website for accurate measurements.

SPLINTER  
FREE EDGES

HIGH TENSILE  
STRENGTH