



Rayon Tipped Applicator- Dry Transport System

Standard double tip, polystyrene handle



Description

Puritan's traditional dry transport tube is perfect for specimen collection and transportation to the lab or use in diagnostic kits. This speciality applicator is made with USP grade rayon fiber that is carded on site to maintain consistency. It is boned securely to a semi-flexible polystyrene handle by an aqueous based adhesive.

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.00" (152.4 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Qualities

- Medical grade quality
- Absorbent tip
- Soft and smooth tip
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Double swab to accomodate multiple sampling sites
- Made in USA

Applications

- Specimen collection for testing and screening
- Cell collection for DNA testing
- General purpose collection swab



Puritan®

Ref #	Packaging
25-806 2PR TT Sterile	24/30/2 (24 boxes per case, 30 tubes per box, 2 swabs per tube)